


**Product Name: Recombinant Human ADAM9 Protein**  
**Catalog #: PEH4012**



## Summary

|                              |                                                                                                                                                                                                               |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                  | ADAM9_HUMAN                                                                                                                                                                                                   |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE                                                                                                                                                           |
| <b>Endotoxin level</b>       | Please contact with the lab for this information                                                                                                                                                              |
| <b>Tag</b>                   | C-His                                                                                                                                                                                                         |
| <b>Accession #</b>           | Q13443                                                                                                                                                                                                        |
| <b>Host</b>                  | E.coli                                                                                                                                                                                                        |
| <b>Species</b>               | Human                                                                                                                                                                                                         |
| <b>Molecular Mass</b>        | 10 kDa                                                                                                                                                                                                        |
| <b>Buffer</b>                | Phosphate buffered saline                                                                                                                                                                                     |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.                                                                                             |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| <b>Sequence</b>              | KRSQTYESDG KNQANPSRQP GSVPRHVSPV TPPREVPIYA NRFAVPTYAA<br>KQPQQFPSRP PPPQPKVSSQ                                                                                                                               |

## image



## Background

|                          |                                                                                                                                                                                                                                    |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Alternative Names</b> | Disintegrin and metalloproteinase domain-containing protein 9 (ADAM 9) (EC 3.4.24.-) (Cellular disintegrin-related protein) (Meltrin-gamma) (Metalloprotease/disintegrin/cysteine-rich protein 9) (Myeloma cell metalloproteinase) |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human ADAM9 Protein**  
**Catalog #: PEH4012**

---

