


Product Name: Recombinant Human MFAP5 Protein
Catalog #: PEH4021



Summary

Name	MFAP5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13361
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GDDVTQATPE TFTEDPNLVN DPATDETVLA VLADIAPSTD DLASLSEKNT TAECWDEKFT CTRLYSVHRP

image



Background

Alternative Names Microfibrillar-associated protein 5 (MFAP-5) (MP25) (Microfibril-associated glycoprotein 2) (MAGP-2)

Note

For Research Use Only , Not for Diagnostic Use