


Product Name: Recombinant Human SIGLEC5 Protein
Catalog #: PEH4026



Summary

Name	SIGL5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	O15389
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YELQVQKSVT VQEGLCVLVP CSFSYPWRSW YSSPLYVYW FRDGEIPYYA EVLATNNPDR RVKPETQGRF RLLGDVQKKN CSLSIGDARM EDTGSYFFRV ERGRDVKYSY QQNKLNLEVT ALIEKPDHF LEPLESGRPT RLSCSLPGSC EAGPPLTFSW TGNALSPLDP ETTRSELT TPRPEDHGTN LTCQMKRQGA

image



Background

Alternative Names	Sialic acid-binding Ig-like lectin 5 (Siglec-5) (CD33 antigen-like 2) (Obesity-binding protein 2) (OB-BP2) (OB-binding protein 2) (CD antigen CD170)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SIGLEC5 Protein
Catalog #: PEH4026

