


Product Name: Recombinant Human NFATC1 Protein
Catalog #: PEH4033



Summary

Name	NFAC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95644
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SPCNKRKYSL NGRQPPYSPH HSPTSPHGS PRVSVTDDSW LGNTTQYTSS AIVAAINALT TDSSLDLGDG VPVKSRTTL

image



Background

Alternative Names Nuclear factor of activated T-cells, cytoplasmic 1 (NF-ATc1) (NFATc1) (NFAT transcription complex cytosolic component) (NF-ATc) (NFATc)

Note

For Research Use Only , Not for Diagnostic Use