


Product Name: Recombinant Human FGF5 Protein
Catalog #: PEH4035



Summary

Name	FGF5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P12034
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YNTYASAIHR TEKTGREWYV ALNKRKGAKR GCSPRVKPQH ISTHFLPRFK QSEQPELSFT VTVPEKKKPP SPIKPKIPLS APRKNTNSVK

image



Background

Alternative Names	Fibroblast growth factor 5 (FGF-5) (Heparin-binding growth factor 5) (HBGF-5) (Smag-82)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use