


Product Name: Recombinant Human ERBB3 Protein
Catalog #: PEH4041



Summary

Name	ERBB3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P21860
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NLGESCQESA VSGSSERCPR PVSLHPMPRG CLASESSEGH VTGSEAEIQE KVSMCRSRSR SRSPRPRGDS AYHSQRHSL TPVTPLSPPG LEEDVNGYV MPDTHLKGTP SSREGTLSSV GLSSVLGTEE EDEDEEY EYM NRRRRHSPPH

image



Background

Alternative Names	Receptor tyrosine-protein kinase erbB-3 (EC 2.7.10.1) (Proto-oncogene-like protein c-ErbB-3) (Tyrosine kinase-type cell surface receptor HER3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use