


Product Name: Recombinant Human NID2 Protein
Catalog #: PEH4059



Summary

Name	NID2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q14112
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NEDNLDYYDV NEEAEYLPG EPEALNGHS SIDVSFQSKV DTKPLESST LDPHTKEGTS LGEVGGPDLK GQVEPWDERE TRSPAPPEVD RDSLAPSWET PPYPENGSI

image



Background

Alternative Names Nidogen-2 (NID-2) (Osteonidogen)

Note

For Research Use Only , Not for Diagnostic Use