

Product Name: Recombinant Human SELPLG Protein
Catalog #: PEH4062



Summary

Name	SELPL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q14242
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YDFLPETEPP EMLRNSTDIT PLTGPGTPES TTVEPAARRS TGLDAGGAVT ELTTELANMG NLSTDSAAME IQTTQPAATE AQTTPVPVTE AQTTPLAATE AQTTRLTATE AQTTPLAATE AQTTPPAATE AQTTPQTGLE

image



Background

Alternative Names	P-selectin glycoprotein ligand 1 (PSGL-1) (Selectin P ligand) (CD antigen CD162)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use