


Product Name: Recombinant Human ADGRG2 Protein
Catalog #: PEH4071



Summary

Name	AGRG2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q8IZP9
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VAKENVRKQW RRYLCCGKLR LAENSDWSKT ATNGLKKQTV NQGVSSSSNS LQSSSNSTNS TTLLVNNDCS VHASGNGNAS TERNGVSFSV QNGDVCLHDF TGKQHMFNEK EDSCNGKGRM ALRRTSKRGS

image



Background

Alternative Names Adhesion G-protein coupled receptor G2 (G-protein coupled receptor 64) (Human epididymis-specific protein 6) (He6)

Note

For Research Use Only , Not for Diagnostic Use