


**Product Name: Recombinant Human ADGRG2 Protein**  
**Catalog #: PEH4072**



## Summary

|                              |                                                                                                                                                                                                                                                                                 |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                  | AGRG2_HUMAN                                                                                                                                                                                                                                                                     |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE                                                                                                                                                                                                                             |
| <b>Endotoxin level</b>       | Please contact with the lab for this information                                                                                                                                                                                                                                |
| <b>Tag</b>                   | N-SUMO & C-His                                                                                                                                                                                                                                                                  |
| <b>Accession #</b>           | Q8IZP9                                                                                                                                                                                                                                                                          |
| <b>Host</b>                  | E.coli                                                                                                                                                                                                                                                                          |
| <b>Species</b>               | Human                                                                                                                                                                                                                                                                           |
| <b>Molecular Mass</b>        | 36 kDa                                                                                                                                                                                                                                                                          |
| <b>Buffer</b>                | Phosphate buffered saline                                                                                                                                                                                                                                                       |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.                                                                                                                                                              |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.                                                                  |
| <b>Sequence</b>              | DTDNSSLSP   PAKLSVVSFA   PSSNGTPEVE   TTSLNDVTLS   LLPSNETEKT<br>KITIVKTFNA   SGVKPQRNIC   NLSSICNDSA   FFRGEIMFQY   DKESTVPQNN<br>HITNGTLTG   LSLSELKRSE   LNKTLQTLSE   TYFIMCATAE   AQSTLNCTFT<br>IKLNNTMNAC   AVIAALERVK   IRPMEHCCCS   VRIPCPSPE   ELEKLQCDLQ<br>DPIVCLADHP |

## image



## Background

|                          |                                                                                                                    |
|--------------------------|--------------------------------------------------------------------------------------------------------------------|
| <b>Alternative Names</b> | Adhesion G-protein coupled receptor G2 (G-protein coupled receptor 64) (Human epididymis-specific protein 6) (He6) |
|--------------------------|--------------------------------------------------------------------------------------------------------------------|

## Note

**Product Name: Recombinant Human ADGRG2 Protein**  
**Catalog #: PEH4072**



---

For Research Use Only , Not for Diagnostic Use