


Product Name: Recombinant Human ST6GALNAC5 Protein
Catalog #: PEH4079



Summary

Name	SIA7E_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BVH7
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QQQQQQQQA SATGSSQPAA ESSTQQRPGV PAGPRPLDGY LGVADHKPLK MHCRCALVT SSGHLLHSRQ GSQIDQTECV IRMNDAPTRG YGRDVGNRST LRVIAHSSIQ

image



Background

Alternative Names Alpha-N-acetylgalactosaminide alpha-2,6-sialyltransferase 5 (EC 2.4.99.-) (GD1 alpha synthase) (GalNAc alpha-2,6-sialyltransferase V) (ST6GalNAc V) (ST6GalNAcV) (Sialyltransferase 7E) (SIAT7-E)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ST6GALNAC5 Protein
Catalog #: PEH4079

