


Summary

Name	CDHR5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q9HBB8
Host	E.coli
Species	Human
Molecular Mass	57 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SVNKDIFEVE ENTNVTEPLV DIHVPEGQEV TLGALSTPFA FRIQGNQLFL NVTPTYEEKS LLEAQLLCQS GGTLVTQLRV FVSVLDVNDN APEFPFKTKE IRVEEDTKVN STVIPETQLQ AEDRDKDDIL FYTLQEMTAG ASDYFSLVSV NRPALRLDRP LDFYERPNMT FWLLVRDTPG ENVEPSHTAT ATLVLNVPVPA DLRPPWFLPC TFSDBGYVICQ AQYHGAVPTG HILPSPLVLR PGPIYAEDGD RGINQPIIYS IFRGNVNGTF IHPDSGNLT VARSVPSPMT FLLLVKGGQQA DLARYSVTQV TVEAVAAAGS PPRFPQRLYR GTVARGAGAG VVVKDAAAPS QPLRIQAQDP EFDLNSAIT YRITNHSHFR MEGEVVLTIT TLAQAGAFYA

image



Background

Alternative Names Cadherin-related family member 5 (Mu-protocadherin) (Mucin and

Product Name: Recombinant Human CDHR5 Protein
Catalog #: PEH4085



cadherin-like protein) (Mucin-like protocadherin) (MLPCDH)

Note

For Research Use Only , Not for Diagnostic Use