

Product Name: Recombinant Human BAIAP2 Protein
Catalog #: PEH4093



Summary

Name	BAIP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UQB8
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GVTGPDGEDY SPWADRKAAQ PKSLSPQSQ SKLSDSYSNT LPVRKSVTPK NSYATTENKT LPRSSSMAAG

image



Background

Alternative Names	Brain-specific angiogenesis inhibitor 1-associated protein 2 (BAI-associated protein 2) (BAI1-associated protein 2) (Protein BAP2) (Fas ligand-associated factor 3) (FLAF3) (Insulin receptor substrate p53/p58) (IRS-58) (IRSp53/58) (Insulin receptor substrate protein of 53 kDa) (IRSp53) (Insulin receptor substrate p53)
--------------------------	--

Note

Product Name: Recombinant Human BAIAP2 Protein
Catalog #: PEH4093



For Research Use Only , Not for Diagnostic Use