


Product Name: Recombinant Human ROS1 Protein
Catalog #: PEH4103



Summary

Name	ROS1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P08922
Host	E.coli
Species	Human
Molecular Mass	52 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VFDVDLEHKV KYPREVKIHN RNSTIISFSV YPLLSRLYWT EVSNFGYQMF YYSIISHTLH RILQPTATNQ QNKRNQCSCN

image



Background

Alternative Names	Proto-oncogene tyrosine-protein kinase ROS (EC 2.7.10.1) (Proto-oncogene c-Ros) (Proto-oncogene c-Ros-1) (Receptor tyrosine kinase c-ros oncogene 1) (c-Ros receptor tyrosine kinase)
--------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Note

For Research Use Only , Not for Diagnostic Use