


Product Name: Recombinant Human CEBPA Protein
Catalog #: PEH4109



Summary

Name	CEBPA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	P49715
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AKAAVGPTGG GGGGDFDYPG APAGPGGAVM PGGAHGPPPG YGCAAAGYLD GRLEPLYERV GAPALRPLVI KQEPREDEA KQLALAGLFP YQPPPPPPPS HHPHPPPPAH LAAPHLQFQI AHCGQTTMHL QPGHPTPPPT PVPSHPAPA LGAAGLPGPG SALKGLGAAH PDLRASGGSG AGKAKKSVDK NSNEYRVRRE RNNIAVRKSR DKAKQRNVET QQKVLELTSNDRLRKRVEQ LSRELDLRLG

image



Background

Alternative Names CCAAT/enhancer-binding protein alpha (C/EBP alpha)

Note

Product Name: Recombinant Human CEBPA Protein
Catalog #: PEH4109



For Research Use Only , Not for Diagnostic Use