


Product Name: Recombinant Human RELT Protein
Catalog #: PEH4112



Summary

Name	TR19L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q969Z4
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DPGQGTLCRP CPPGTFSAAW GSSPCQPHAR CSLWRRLEAQ VGMATRDTLC GDCWPGWFGP WGVPRVPCQP CSWAPLGTHG CDEWGRRARR GVEVAAGASS GGETRQPGNG TRAGGPEETA

image



Background

Alternative Names Tumor necrosis factor receptor superfamily member 19L (Receptor expressed in lymphoid tissues)

Note

For Research Use Only , Not for Diagnostic Use