


Product Name: Recombinant Human ADAMTS1 Protein
Catalog #: PEH4131



Summary

Name	ATS1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UHI8
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TPCSPDSTSV CVQGQCVKAG CDRIIDSKKK FDKCGVCGGN GSTCKKISGS VTSAKPGYHD IITPTGATN IEVKQRNQRG

image



Background

Alternative Names A disintegrin and metalloproteinase with thrombospondin motifs 1 (ADAM-TS 1) (ADAM-TS1) (ADAMTS-1) (EC 3.4.24.-) (METH-1)

Note

For Research Use Only , Not for Diagnostic Use