


Product Name: Recombinant Human ST6GAL1 Protein
Catalog #: PEH4134



Summary

Name	SIAT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P15907
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KGSYYDSFKL QTKEFQVLKS LGKLAMGSDS QSVSSSTQD PHRGRQTLGS LRGLAKAKPE ASFQVWNKDS

image



Background

Alternative Names Beta-galactoside alpha-2,6-sialyltransferase 1 (Alpha 2,6-ST 1) (EC 2.4.99.1) (B-cell antigen CD75) (CMP-N-acetylneuraminase-beta-galactosamide-alpha-2,6-sialyltransferase 1) (ST6Gal I) (ST6GalII) (Sialyltransferase 1)

Note

For Research Use Only , Not for Diagnostic Use