


Product Name: Recombinant Human VCPKMT Protein
Catalog #: PEH4140



Summary

Name	MT21D_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H867
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DVVVTDLEEL QDLLKMNINM NKHLVTGSVQ AKVLKWGEEI EGFSPPDFI LMADCIYEE SLEPLLKTLK DISGFETCII CCYEQRTMGK NPEIEKKYFE LLQLDFDFEK IPLEKHDEEY

image



Background

Alternative Names	Protein-lysine methyltransferase METTL21D (EC 2.1.1.-) (Methyltransferase-like protein 21D) (VCP lysine methyltransferase) (VCP-KMT) (Valosin-containing protein lysine methyltransferase)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human VCPKMT Protein
Catalog #: PEH4140

