


Product Name: Recombinant Human IL12RB2 Protein
Catalog #: PEH4153



Summary

Name	I12R2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99665
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DPEPLVISEV LHQVTPVFRH PPCSNWPQRE KGIQGHQASE KDMMHSAASP PPPRALQAES RQLVDLYKVL ESRGSDPKPE

image



Background

Alternative Names	Interleukin-12 receptor subunit beta-2 (IL-12 receptor subunit beta-2) (IL-12R subunit beta-2) (IL-12R-beta-2) (IL-12RB2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use