


Product Name: Recombinant Human NFKBIE Protein
Catalog #: PEH4154



Summary

Name	IKBE_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O00221
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PPRAGRPLCQ ADPQPGKAAR RSLEPDPAQT GPRPARAAGM SEARKGPDEA EESQYDSGIE SLRSLRSLPE

image



Background

Alternative Names	NF-kappa-B inhibitor epsilon (NF-kappa-BIE) (I-kappa-B-epsilon) (IkB-E) (IkB-epsilon) (IkappaBepsilon)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use