


Product Name: Recombinant Human LRIG3 Protein
Catalog #: PEH4156



Summary

Name	LRIG3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6UXM1
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SITNTDETNL PADIPSYLSS QGTLADRQDG YVSESSESHH QFVTSSGAGF FLPQHDSST CHIDNSSEAD VEAATDLFLC PFLGSTGPMY LKGNVYGSDP FETYHTGCSP DPRTVLMHDY EPSYIKKKEC YPCSHPEES CERFSNISW PSHVRKLLNT SYSHNEGPGM KNLCLNKSSL DFSANPEPAS VASSNSFMGT FGKALRRPHL DAYSSFGQPS DCQPRAFYLK AHSSPDLDG SEEDGKERTD FQEENHICTF KQTLENYRTP

image



Background

Alternative Names Leucine-rich repeats and immunoglobulin-like domains protein 3 (LIG-3)

Product Name: Recombinant Human LRIG3 Protein
Catalog #: PEH4156



Note

For Research Use Only , Not for Diagnostic Use