


Product Name: Recombinant Human PNPT1 Protein
Catalog #: PEH4157



Summary

Name	PNPT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TCS8
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KTISKPRASR KENGPVVETV QVPLSKRAKF VGPGGYNLKK LQAETGVTIS QVDEETFSVF APTPSAMHEA RDFITEICKD

image



Background

Alternative Names	Polyribonucleotide nucleotidyltransferase 1, mitochondrial (EC 2.7.7.8) (3'-5' RNA exonuclease OLD35) (PNPase old-35) (Polynucleotide phosphorylase 1) (PNPase 1) (Polynucleotide phosphorylase-like protein)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use