


Product Name: Recombinant Human FABP5 Protein
Catalog #: PEH4161



Summary

Name	FABP5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q01469
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ATVQQLLEGRW RLVDSKGFDE YMKELGVGIA LRKMGAMAKP DCIITCDGKN LTIKTESTLK TTQFSCTLGE KFEETTADGR KTQTVCNFTD GALVQHQEWD GKESTITRKL KDGKLVVECV MNNVTCTRIY EKVE

image



Background

Alternative Names	Fatty acid-binding protein 5 (Epidermal-type fatty acid-binding protein) (E-FABP) (Fatty acid-binding protein, epidermal) (Psoriasis-associated fatty acid-binding protein homolog) (PA-FABP)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FABP5 Protein
Catalog #: PEH4161

