


Product Name: Recombinant Human FXN Protein
Catalog #: PEH4174



Summary

Name	FRDA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q16595
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ETLDSLAEFF EDLADKPYTF EDYDVSGSG VLTVKLGGDL GTYVINKQTP NKQIWLSSPS SGPKRYDWTG

image



Background

Alternative Names Frataxin, mitochondrial (EC 1.16.3.1) (Friedreich ataxia protein) (Fxn) [Cleaved into: Frataxin intermediate form (i-FXN); Frataxin(56-210) (m56-FXN); Frataxin(78-210) (d-FXN) (m78-FXN); Frataxin mature form (Frataxin(81-210)) (m81-FXN)]

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FXN Protein
Catalog #: PEH4174

