


**Product Name: Recombinant Human MMP15 Protein**  
**Catalog #: PEH4176**



## Summary

<b>Name</b>	MMP15_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P51511
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	10 kDa
<b>Buffer</b>	Lyophilized from Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	LPEDDLRGIQ QLYGTPDGQP QPTQPLPTVT PRRPGRPDHR PRRPPQPPPP GGKPERPPKP GPPVQPRATE

## image



## Background

**Alternative Names** Matrix metalloproteinase-15 (MMP-15) (EC 3.4.24.-) (Membrane-type matrix metalloproteinase 2) (MT-MMP 2) (MTMMP2) (Membrane-type-2 matrix metalloproteinase) (MT2-MMP) (MT2MMP) (SMCP-2)

## Note

For Research Use Only , Not for Diagnostic Use