


Product Name: Recombinant Human SCARA5 Protein
Catalog #: PEH4189



Summary

Name	SCAR5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q6ZMJ2
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KGEKGEKGD R AGDASGVEAP MMIRLVNGSG PHEGRVEVYH DRRWGTVCDD GWDKKGDDVV CRMLGFRGVE EVYRTARFGQ GTGRIWMDDV ACKGTEETIF RCSFSKWGVT NCGHAEDASV

image



Background

Alternative Names Scavenger receptor class A member 5 (Scavenger receptor hlg)

Note

For Research Use Only , Not for Diagnostic Use