


Product Name: Recombinant Human XDH Protein
Catalog #: PEH4191



Summary

Name	XDH_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P47989
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NLEDKCGKLD PTFASATLLF QKDPPADVQL FQEVPKGQSE EDMVGRPLPH LAADMQASGE AVYCDDIPRY ENELSLRLVT STRAHAKIKS

image



Background

Alternative Names Xanthine dehydrogenase/oxidase [Includes: Xanthine dehydrogenase (XD) (EC 1.17.1.4); Xanthine oxidase (XO) (EC 1.17.3.2) (Xanthine oxidoreductase) (XOR)]

Note

For Research Use Only , Not for Diagnostic Use