


**Product Name: Recombinant Human PGR Protein**  
**Catalog #: PEH4199**



## Summary

<b>Name</b>	PRGR_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P06401
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	20 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	IPISLDGLLF PRPCQGQDPS DEKTQDQQL SDVEGAYSRA EATRGAGGSS SSPPEKDSGL LDSVLDTLA PSGPGQSQPS PPACEVTSSW CLFGPELPED PPAAPATQRV LSPLMSRSGC KVGDSSTAA AHKVLPRGLS PARQLLLPAS ESPHWSGAPV KPSPQAAAVE

## image



## Background

**Alternative Names** Progesterone receptor (PR) (Nuclear receptor subfamily 3 group C member 3)

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human PGR Protein**  
**Catalog #: PEH4199**

---

