


Product Name: Recombinant Human PRDX5 Protein
Catalog #: PEH4207



Summary

Name	PRDX5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P30044
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAPIKVGDAI PAVEVFEGEP GNKVNLAELF KGKKGVLFV PGAFTPGCSK THLPGFVEQA EALKAKGVQV VACLSVNDAF VTGEWGRAHK AEGKVRLAD PTGAFGKETD LLLDDSLVSI FGNRRRLKRF MVVQDGIVKA LNVEPDGTGL TCSLAPNIIS QL

image



Background

Alternative Names	Peroxiredoxin-5, mitochondrial (EC 1.11.1.24) (Alu corepressor 1) (Antioxidant enzyme B166) (AOEB166) (Liver tissue 2D-page spot 71B) (PLP) (Peroxiredoxin V) (Prx-V) (Peroxisomal antioxidant enzyme) (TPx type VI) (Thioredoxin peroxidase PMP20) (Thioredoxin-dependent peroxiredoxin 5)
--------------------------	---

Product Name: Recombinant Human PRDX5 Protein
Catalog #: PEH4207



Note

For Research Use Only , Not for Diagnostic Use