


Product Name: Recombinant Human FHL1 Protein
Catalog #: PEH4215



Summary

Name	FHL1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13642
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HCVKCNKAIT SGGITYQDQP WHADCFVCVT CSKLAGQRF TAVEDQYYCV DCYKNFVAKK CAGCKNPITG KRTVSRVSHV VSKARKPPVC

image



Background

Alternative Names Four and a half LIM domains protein 1 (FHL-1) (Skeletal muscle LIM-protein 1) (SLIM) (SLIM-1)

Note

For Research Use Only , Not for Diagnostic Use