


**Product Name: Recombinant Human DECR1 Protein**  
**Catalog #: PEH4217**



## Summary

<b>Name</b>	DECR_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q16698
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	12 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	GPIKTKGAFS RLDPTGTFEK EMIGRIPCGR LGTVEELANL AAFLCSDYAS WINGAVIKFD GGEEVLISGE FNDLRKVTKE QWDTIEELIR

## image



## Background

<b>Alternative Names</b>	2,4-dienoyl-CoA reductase, mitochondrial (EC 1.3.1.34) (2,4-dienoyl-CoA reductase [NADPH]) (4-enoyl-CoA reductase [NADPH]) (Short chain dehydrogenase/reductase family 18C member 1)
--------------------------	--

## Note

For Research Use Only , Not for Diagnostic Use