



Product Name: Recombinant Human MEGF9 Protein
Catalog #: PEH4222



Summary

Name	MEGF9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9H1U4
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VTGGGGAAGQ VDASPGGLR GEP SHPFRA TAPTAQAPRT GPPRATVHRP LAATSPAQSP ETTPLWATAG PSSTTFQAPL GPSPTTPPAA ERTSTTSQAP TRPAP TTLST TTGPAPTTPV ATTVPAPTTP RTPTDLPSS SNSSVLPTPP ATEAPSSPPP EYVCNCSVVG SLNVNRCNQT TGQCECRPGY QGLHCETCKE GFYLNYSGL

image



Background

Alternative Names	Multiple epidermal growth factor-like domains protein 9 (Multiple EGF-like domains protein 9) (Epidermal growth factor-like protein 5) (EGF-like protein 5)
--------------------------	---

Product Name: Recombinant Human MEGF9 Protein
Catalog #: PEH4222



Note

For Research Use Only , Not for Diagnostic Use