


Product Name: Recombinant Human STAB1 Protein
Catalog #: PEH4225



Summary

Name	STAB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NY15
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TPGDPKRTIG QILASTEAFS RFETILENCG LPSILDGPGP FTVFAPSNEA VDSLRDGRLI YLFTAGLSKL QELVRYHIYN HGQLTVEKLI SKGRILTMAN QVLAVNISEE GRILLGPEGV PLQRVDVMAA NGVIHMLDGI LLPPTILPIL PKHCSEEQHK IVAGSCVDCQ ALNTSTCPPN SVKLDIFPKE CVYIHDPTGL NVLKKGCASY

image



Background

Alternative Names	Stabilin-1 (Fasciclin, EGF-like, laminin-type EGF-like and link domain-containing scavenger receptor 1) (FEEL-1) (MS-1 antigen)
--------------------------	---

Note

Product Name: Recombinant Human STAB1 Protein
Catalog #: PEH4225



For Research Use Only , Not for Diagnostic Use