


## Summary

<b>Name</b>	NUDT5_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9UKK9
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	25 kDa
<b>Buffer</b>	Lyophilized from Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	MESQEPTES QNGKQYIISE ELISEGKWVK LEKTTYMDPT GKTRTWESVK RTRRKEQTAD GVAVIPVLQR TLHYECIVLV KQFRPPMGY CIEFPAGLID DGETPEAAAL RELEEETGYK GDIAECSPAV CMDPGLSNT IHIVTVTING DDAENARPKP KPGDGEFVEV ISLPKNDLLQ RLDALVAEEH LTVDARVYSY ALALKHANAK PFEVPFLKF

## image



## Background

<b>Alternative Names</b>	ADP-sugar pyrophosphatase (EC 3.6.1.13) (8-oxo-dGDP phosphatase) (EC 3.6.1.58) (Nuclear ATP-synthesis protein NUDIX5) (EC 2.7.7.96) (Nucleoside diphosphate-linked moiety X motif 5) (Nudix motif 5) (hNUDT5) (YSA1H)
--------------------------	---

**Product Name: Recombinant Human NUDT5 Protein**  
**Catalog #: PEH4226**



---

**Note**

For Research Use Only , Not for Diagnostic Use