


Product Name: Recombinant Human CORO1A Protein
Catalog #: PEH4234



Summary

Name	COR1A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P31146
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GGRDAGPLLI SLKDGYPVK SRELVRNRL DTGRRRAPE ASGTPSSDAV SRLEEEMRKL QATVQELQKR LDRLEETVQA

image



Background

Alternative Names Coronin-1A (Coronin-like protein A) (Clipin-A) (Coronin-like protein p57) (Tryptophan aspartate-containing coat protein) (TACO)

Note

For Research Use Only , Not for Diagnostic Use