


Product Name: Recombinant Human ZBP1 Protein
Catalog #: PEH4235



Summary

Name	ZBP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H171
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SISPGVAGPG GVAGSGEGEP GEDAGRRPAD TQSRSHFPRD IGQPITPSHS KLTPKLETMT LGNRSHKAAE GSHYVDEASH

image



Background

Alternative Names Z-DNA-binding protein 1 (DNA-dependent activator of IFN-regulatory factors) (DAI) (Tumor stroma and activated macrophage protein DLM-1)

Note

For Research Use Only , Not for Diagnostic Use