


Product Name: Recombinant Human KCNIP4 Protein
Catalog #: PEH4238



Summary

Name	KCIP4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6PIL6
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LNWAFNLYDI NKDGYITKEE MLDIMKAIYD MMGKCTYPVL KEDAPRQHVE TFFQKMDKNK DGVVTIDEFI

image



Background

Alternative Names	Kv channel-interacting protein 4 (KChIP4) (A-type potassium channel modulatory protein 4) (Calsenilin-like protein) (Potassium channel-interacting protein 4)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use