



Product Name: Recombinant Human TYMP Protein
Catalog #: PEH4241



Summary

Name	TYPH_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P19971
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VGDKVSLVLA PALAACGCKV PMISGRGLGH TGGTLDKLES IPGFNVIQSP EQMQVLLDQA GCCIVGQSEQ

image



Background

Alternative Names	Thymidine phosphorylase (TP) (EC 2.4.2.4) (Gliostatin) (Platelet-derived endothelial cell growth factor) (PD-ECGF) (TdRPase)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use