


Product Name: Recombinant Human HS3ST3B1 Protein
Catalog #: PEH4248



Summary

Name	HS3SB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9Y662
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PDGTPPRLPF RAPPATPLAS GKEMAEGAAS PEEQSPEVPD SPSPISFFS GSGSKQLPQA IIIGVKKGGT RALLEFLRVH PDVRAVGAEP HFFDRSYDKG

image



Background

Alternative Names Heparan sulfate glucosamine 3-O-sulfotransferase 3B1 (EC 2.8.2.30)
(Heparan sulfate D-glucosaminyl 3-O-sulfotransferase 3B1) (3-OST-3B)
(Heparan sulfate 3-O-sulfotransferase 3B1) (h3-OST-3B)

Note

For Research Use Only , Not for Diagnostic Use