


Product Name: Recombinant Human CEACAM8 Protein
Catalog #: PEH4249



Summary

Name	CEAM8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P31997
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SNNSNPVEDK DAVAFTCEPE TQNTTYLWWV NGQSLPVSPR LQLSNGNRTL TLLSVTRNDV GPYECEIQNP ASANFSDPVT

image



Background

Alternative Names	Carcinoembryonic antigen-related cell adhesion molecule 8 (CD67 antigen) (Carcinoembryonic antigen CGM6) (Non-specific cross-reacting antigen NCA-95) (CD antigen CD66b)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use