


Product Name: Recombinant Human BCAM Protein
Catalog #: PEH4268



Summary

Name	BCAM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P50895
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VGSPSTPAGW VREGDTVQLL CRGDGSPSPE YTLFRLQDEQ EEVLNVNLEG NLTLEGVTRG QSGTYGCRVE DYDAADDVQL SKTLELRVAY

image



Background

Alternative Names	Basal cell adhesion molecule (Auburger B antigen) (B-CAM cell surface glycoprotein) (F8/G253 antigen) (Lutheran antigen) (Lutheran blood group glycoprotein) (CD antigen CD239)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use