


Product Name: Recombinant Human FGF16 Protein
Catalog #: PEH4270



Summary

Name	FGF16_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	O43320
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LDWDLHGFSS SLGNVPLADS PGFLNERLGQ IEGKLQRGSP TDFAHKLGIL RRRQLYCRTG FHLEIFPNGT VHGTRHDHSR FGILEFISLA VGLISIRGVD SGLYLGMNER GELYGSKKLT RECVFREQFE ENWYNTYAST LYKHSDSERQ YYVALNKDGS PREGYRTKRH QKFTHFLPRP VDPSKLPSMS

image



Background

Alternative Names Fibroblast growth factor 16 (FGF-16)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FGF16 Protein
Catalog #: PEH4270

