


Product Name: Recombinant Human IL24 Protein
Catalog #: PEH4277



Summary

Name	IL24_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q13007
Host	E.coli
Species	Human
Molecular Mass	52 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SDAESCYLEVH TLLFYLKTV FKNYHNRTVE VRTLKSFSTL ANNFVLIVSQ LQPSQENEMF SIRDSEHRRF LLFRRAFKQL

image



Background

Alternative Names	Interleukin-24 (IL-24) (Melanoma differentiation-associated gene 7 protein) (MDA-7) (Suppression of tumorigenicity 16 protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use