


Product Name: Recombinant Human CEACAM3 Protein
Catalog #: PEH4280



Summary

Name	CEAM3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P40198
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MPLSVAEGKE VLLLVHNL PQ HLF GYSWYKG ERVDGNSLIV GYVIGTQQAT PGAAYSGRET IYTNASLLIQ NVTQNDIGFY TLQVIKSDLV NEEATGQFHV YQENAPGLPV

image



Background

Alternative Names	Carcinoembryonic antigen-related cell adhesion molecule 3 (Carcinoembryonic antigen CGM1) (CD antigen CD66d)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use