


Product Name: Recombinant Human EMG1 Protein
Catalog #: PEH4287



Summary

Name	NEP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92979
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RERSGGEEQAQ DWDALPPKRP RLGAGNKIGG RRLIVVLEGA SLETVKVGKT YELLNCDKHK SILLKNGRDP GEARPDITHQ

image



Background

Alternative Names	Ribosomal RNA small subunit methyltransferase NEP1 (EC 2.1.1.-) (18S rRNA (pseudouridine(1248)-N1)-methyltransferase) (18S rRNA Psi1248 methyltransferase) (Nucleolar protein EMG1 homolog) (Protein C2f) (Ribosome biogenesis protein NEP1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human EMG1 Protein
Catalog #: PEH4287

