


Product Name: Recombinant Human COQ6 Protein
Catalog #: PEH4290



Summary

Name	COQ6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y2Z9
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AVVATLHLS E ATENNAWQR FLPSGPIALL PLSDTLSSLV WSTSHEHAAE LVSMDEEK FV DAVNSAFWSD ADHTDFIDTA

image



Background

Alternative Names Ubiquinone biosynthesis monooxygenase COQ6, mitochondrial (EC 1.14.13.-) (Coenzyme Q10 monooxygenase 6)

Note

For Research Use Only , Not for Diagnostic Use