


Product Name: Recombinant Human RARRES1 Protein
Catalog #: PEH4303



Summary

| | |
|------------------------------|--|
| Name | TIG1_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | N-SUMO & N-His |
| Accession # | P49788 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 41 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | DAGVPRLLQ QAARAALHFF NFRSGSPSAL RVLAEVQEGR AWINPKEGCK VHVVFSTERY NPESLLQEGE GRLGKCSARV FFKNQKPRPT INVCTRLIE KKKRQQEDYL LYKQMKQLKN PLEIVSIPDN HGHIDPSLRL IWDLAFLGSS YVMWEMTTQV SHYYLAQLTS VRQWKTNDDT IDFDYTVLLH ELSTQEIIPC RIHLVWYPGK PLKVKYHCQE LQTPEEASGT EEGSAVVPT |

image



Background

| | |
|--------------------------|--|
| Alternative Names | Retinoic acid receptor responder protein 1 (Phorbol ester-induced gene 1 protein) (PERG-1) (RAR-responsive protein TIG1) (Tazarotene-induced gene 1 protein) |
|--------------------------|--|

Product Name: Recombinant Human RARRES1 Protein
Catalog #: PEH4303



Note

For Research Use Only , Not for Diagnostic Use