


Product Name: Recombinant Human CLEC10A Protein
Catalog #: PEH4308



Summary

Name	CLC10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IUN9
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TLRTDFSNT SNTVAEIQAL TSQGSSLEET IASLKAIEVEG FKQERQAGVS ELQEHTTQKA HLGHCPCPS

image



Background

Alternative Names C-type lectin domain family 10 member A (C-type lectin superfamily member 14) (Macrophage lectin 2) (CD antigen CD301)

Note

For Research Use Only , Not for Diagnostic Use