


Product Name: Recombinant Human PON2 Protein
Catalog #: PEH4311



Summary

Name	PON2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15165
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SREVESVDLP HCHLIKGIEA GSEDIDILPN GLAFFSVGLK FPGLHSFAPD KPGGILMMDL KEEKPRAREL RISRGFDLAS FNPHGISTFI DNDDTVYLFV VNHPEFKNTV EIFKFEEAEN SLLHLKTVKH ELLPSVNDIT AVGPAHFYAT NDHYFSDPFL KYLETYLN LH WANVVYSPN EVKVVAEGFD SANGINISPD DKYIYVADIL AHEIHVLEKH TNMNLTLQKV LELDTLVDNL SIDPSSGDIW VGCHPNGQKL FVYDPNNPPS SEVLRIQNIL SEKPTVTTVY ANNGSVLQGS SVASVYDGKL LIGTLYHRAL

image



Background

Alternative Names Serum paraoxonase/arylesterase 2 (PON 2) (EC 3.1.1.2) (EC 3.1.1.81) (Aromatic esterase 2) (A-esterase 2) (Serum aryldialkylphosphatase 2)

Product Name: Recombinant Human PON2 Protein
Catalog #: PEH4311



Note

For Research Use Only , Not for Diagnostic Use